



AIVIA STAND

For Automated
External
Defibrillators
(AED)



AIVIA[®] STAND

OVERVIEW

The AIVIA Stand was designed for the AIVIA cabinet brand and utilized in public access areas in applications when there is no structure available to install cabinet to.


AIVIA Stand Features:

- Heavy duty aluminum and highly visible red color
- Suitable for indoor/outdoor use in any application or location
- Stand provides a structure to attach cabinet to in open access areas
- Locked storage space
- Protective case inside locked storage space to hold power supply
- 2 heavy duty keys
- (4) Screws to attach cabinet
- Holes in bottom to secure to ground
- Electrical wiring can be wired from cabinet and through the stand either in back or at bottom



LABELING

- Visible AED signage



DIMENSIONS & WEIGHT

- Width of stand : 15 9/32" (1)
- Space on each side of stand for application clearance: 4 11/32" (2)
- Space above stand for application clearance: 3 15/16"
- Height from top of stand to bottom: 68 5/16 "
- Height from bottom of stand to bottom of cabinet: 47 1/4"
- Depth of stand: 7 13/16"
- Depth of stand with door ajar: 55 1/8"
- Weight: 30lbs

